


KANEKA REN



Balloon Dilation Catheter

- Radial Balloon Dilation
- Seamless Tip Design
- Visual Marker





Radial Balloon Dilation


Radial Balloon Dilation prevents slipping from vater papilla and provide withdrawal resistance

Visual Marker

Visual marker design for easier positioning

Seamless Tip Design


Smooth insertion of the catheter into narrow parts



Radiopaque Marker

Radiopaque marker for easier fluoroscopic guidance of balloon positioning

Specifications



Ordering Information

Product Code	Balloon Diameter	Balloon Length	Accepted GW	Sheath Length	Max Pressure	Working Channel
RN25-0320-18	3mm	20mm	0.025"	1,800mm	8atm	2.8mm
RN25-0430-18	4mm	30mm				
RN25-0630-18	6mm					
RN25-0830-18	8mm	40mm				
RN25-1240-18	10-11-12mm					
RN25-1540-18	13-14-15mm					
RN25-1840-18	16-17-18mm					

*Please refer to IFU before use